


Weitere Informationen siehe unter [www.ostene.com](http://www.ostene.com)

		Avitene Instat Ultrafoam	Gelfoam Spongostan Surgifoam	Thrombin	FloSeal	Surgicel	Bone Wax
Agent	Alkylene Oxide Copolymers	Microfibrillar Collagen	Gelatin	Thrombin	Gelatin granules Thrombin	Oxidized Regenerated Cellulose	Beeswax
Supplied	Bars	Sheets/powder	Sheets/powder	Powder / Solution	Gel in Syringe	Sheets	Bars
Ingredients	Water soluble Synthetic copolymers	Bovine Corium Collagen	Animal Gelatin	Bovine Thrombin	Bovine Gelatin Bovine Thrombin	Wood pulp	Beeswax, Vaseline, Isopropylpalmitate
Action	Mechanical occlusion	Platelet release stimulation	Mechanical occlusion, swells, Gel formation	Activates Fibrin	Mechanical + Activates Fibrin	Mechanical occlusion Swells, gel formation	Mechanical occlusion
Remains at Site	24 - 48 hours	<12 weeks	4 - 6 weeks	minutes	6 - 8 weeks	Variable	Indefinitely
Advantages	<ul style="list-style-type: none"> <li>• Immediate hemostasis</li> <li>• Ease of use</li> <li>• Adheres to wet surfaces</li> </ul>	<p><b>For more details visit: <a href="http://www.ostene.com">www.ostene.com</a></b></p>					
Adverse Effects	<ul style="list-style-type: none"> <li>• None known</li> </ul>						